

# Bioss HRP Polymer Ready-To-Use IHC Kits

## Simplify Staining with Enhanced Sensitivity and Efficient Amplification

**Bioss HRP Polymer Ready-To-Use IHC Kits** feature a HRP polymer-conjugated secondary antibody, which significantly addresses limitations commonly experienced with traditional immunohistochemistry methods. The HRP polymer amplification method **enhances signal detection**, while reducing background noise allowing for high-quality staining even with low-abundance antigens or suboptimal antibody-antigen interactions. With **less primary antibody needed** and a **shortened secondary antibody incubation period**, Bioss HRP Polymer Ready-To-Use IHC Kits enable a staining process that is simpler, faster and more efficient.

## Kit Components

Each kit includes the following **ready-to-use reagents**:

- ✓ HRP Polymer Conjugated Secondary Antibody
- ✓ Endogenous Peroxidase Blocking Solution
- ✓ Normal Goat Serum

### TARGET

### CATALOG NO.

Goat Anti-Mouse & Rabbit IgG H&L (HRP polymer) Ready-to-use Kit	PV-1022
Goat Anti-Rabbit IgG H&L (HRP polymer) Ready-to-use Kit	PV-1023
Goat Anti-Mouse IgG H&L (HRP polymer) Ready-to-use Kit	PV-1024

Each kit is offered in 6ml and 50ml sizes.

To view our full catalog, please visit [Biossusa.com](http://Biossusa.com).